



## Assay Methods Using Peptidyl-AMC (MCA) Substrates (1)


### Principle



The initial rate of increase in the AMC concentration can be monitored 1) fluorometrically at  $\lambda_{\text{ex}} = 380 \text{ nm}$  and  $\lambda_{\text{em}} = 460 \text{ nm}$  (Fig. 1a) or 2) photometrically at  $370 \text{ nm}$  (Fig. 1b).



**Fig. 1a Fluorescence Spectra**



**Fig. 1b UV-Absorption Spectra**


### Reagents

1. Substrate stock solution: Vial, in DMSO at 10 mM
2. AMC stock solution: Content of vial (Code 3099-vAMC), in DMSO at 1 mM
3. Buffer
4. Enzyme solution

### Procedure

Choose the proper conditions for the measurement, such as substrate concentration and sensitivity setting, depending on the purpose of the experiment and the instrument available. Described here is one of the recommended procedures for the fluorometric method (initial-rate method).

1. Set a fluorescence spectrophotometer at  $\lambda_{ex} = 380 \text{ nm}$  and  $\lambda_{em} = 460 \text{ nm}$  at  $25 \text{ }^\circ\text{C}$  (1.0 Relative fluorescence unit at  $10^{-6} \text{ M}$  of AMC)
2. Pipette 2940  $\mu\text{l}$  of buffer and 30  $\mu\text{l}$  of substrate stock solution into the cuvette
3. Incubate in the fluorescence spectrophotometer for 3-4 min (for temperature equilibration)
4. Add 30  $\mu\text{l}$  of enzyme solution
5. Record the increase of the fluorescence intensity for 3-4 min
6. Calculate the amount of AMC released using the following equation



Amount (nmol) of AMC released/min

$$\begin{aligned}
 &= \frac{1 \text{ nmol} \times dF\% \times 3 \text{ ml}}{1 \text{ ml} \times 100\% \times 1 \text{ min}} \\
 &= 0.03 \times dF \text{ nmol/min}
 \end{aligned}$$

\* Photometric measurement can be carried out by the same procedure as that of the fluorometric method using a UV spectrophotometer. Set the wavelength at 370 nm ( $\epsilon$  7700).