

Assay Methods Using Peptidyl-AMC (MCA) Substrates (3)

(Example for Caspase-7 using an Auto-Fluorescence Spectrophotometer for Multiplate)


Reagents

1. Substrate stock solution: Content of vial (Code 3171-v Ac-DEVD-AMC), in DMSO at 10 mM
2. AMC stock solution: Content of vial (Code 3099-vAMC), in DMSO at 1 mM
3. Buffer: ICE standard buffer (100 mM HEPES-KOH, pH 7.5, 10% sucrose (w/v), 0.1% CHAPS (w/v), 10 mM DTT, 0.1 mg/ml ovalbumin)
4. Enzyme solution: rec. Caspase-7 in buffer

Procedure

- a) Initial rate assay for 30 min at $\lambda_{ex} = 390$ nm, $\lambda_{em} = 460$ nm
- b) End-point (30 min) assay at 5 different substrate concentrations at $\lambda_{ex} = 355$ nm or 380 nm, $\lambda_{em} = 460$ nm

1. Set a fluorescence spectropotometer at $\lambda_{ex} = 380$ nm (or 380 nm, 355 nm) and $\lambda_{em} = 460$ nm
Relative fluorescence is determined with 10^{-6} M of AMC at each condition
2. Substrate: Dilute the stock solution with the buffer (x 10, x 20, x 40, x 80, x 160)
3. Caspase-7: Dissolve in buffer
4. Pipette 160 μ l of buffer and 20 μ l of substrate solutions (1, $\frac{1}{2}$, $\frac{1}{4}$, $\frac{1}{8}$, $\frac{1}{16}$ mM) in each well
5. Incubate in the fluorescence spectrophotometer for 3-4 min (for temperature equilibration)
6. Take the multiplate out and add 20 μ l of enzyme solution in each well
7. Mount the plate in the fluorescence spectrometer
- 8 Calculate the amount of released AMC



Results

